



SARS-CoV-2 (N) Ig ELISA Kit (SCSE-Ig)

Catalog Number: EL2045-01

I. Application

This immunoassay kit allows for the qualitative determination of 2019 nCoV-Ig antibody in serum, plasma, saliva and nasal fluid.

- **Size:** 96T
- **Principle:** Sandwich
- **Reactivity:** Universal
- **Species:** Human and Animal

2. Storage and Kit Contents

Storage: 2-8°C for 6 months.


Kit Components	Specifications (96T)	Storage
Coated assay plate	1 plate	2-8°C
Negative Control (Ready-to-use)	1 vial	2-8°C
Positive Control (Ready-to-use)	1 vial	2-8°C
Sample Dilution Buffer	1 bottle	2-8°C
Biotin-conjugated Nucleocapsid (Concentrated)	1 vial	2-8°C(avoid direct light)
Antigen Dilution Buffer	1 bottle	2-8°C
HRP-Streptavidin Conjugate (SABC)	1 vial	2-8°C(avoid direct light)
SABC Dilution Buffer	1 bottle	2-8°C
Wash Buffer (25 x concentrate)	1 bottle	2-8°C
TMB Substrate	1 bottle	2-8°C(avoid direct light)
Stop solution	1 bottle	2-8°C
Plate Sealer	5 pieces	
Instruction manual	1 copy	

NOTE: FOR RESEARCH USE ONLY.

3. Operation Procedure

Principle of the Assay

This kit was based on sandwich enzyme-linked immune-sorbent assay technology. Recombinant 2019-nCoV Nucleocapsid protein (antigen) was pre-coated onto 96-well plates. The Controls, test samples and Biotin-labeled antigen were added to the wells subsequently, and wash with wash buffer. HRP-Streptavidin Conjugate is added and unbound conjugates were washed away with wash buffer. TMB was catalyzed by HRP to produce a blue color product that changed into yellow after adding acidic stop solution. The optical density of developed color is read with a suitable photometer at 450nm with a selected reference wavelength within 650 nm.



4. Sequence of Nucleocapsid (antigen)

MSDNGPQNQRNAPRITFGGSPDSTGSNQNGERSGARSKQRRPQGLPNNTASWFTALT
QHGKEDLKFPRGQGVPIINTSSPDDQIGYYRRATRRIRGGDGKMKDLSRWYFYLLGTG
PEAGLPYGANKDGIIWVATEGALNTPKDHIGTRNPANNAIIVLQLPQGTTLPKGFYAEGS
RGGSQASSRSSSRNRNSTRNTPGSSRGTSPARMAGNGGDAALALLLDRLNQLKESKMSG
KGQQQQGQTVTKKSAAEASKKPRQKRTATKAYNVTQAFGRRGPEQTQGNFGDQELIRQ
GTDYKHWPQIAQFAPSASAFFGMSRIGMEVTPSGTWLTYTGAIKLDDKDPNFKDQVILLN
KHIDAYKTFPPTPEPKDKKKKADETQALPQRQKKQQTVTLPAADLDDFSKQLQQSMSSA
DSTQA

5. Precautions

1. To inspect the validity of experiment operation and the appropriateness of sample dilution proportion, pilot experiment using standards and a small number of samples is recommended.
2. After opening and before using, keep plate dry.
3. Before using the kit, spin tubes and bring down all components to the bottom of tubes.
4. Storage TMB reagents avoid light.
5. Washing process is very important, not fully washing can cause a false positive and high background.
6. Duplicate well assays are recommended for both standard and sample testing.
7. Don't let microplate dry as a dry plate will inactivate the active components on plate.
8. Don't reuse tips and tubes to avoid cross contamination.
9. Avoid using the reagents from different batches together.

6. User Supplied Materials

- Microplate reader (wavelength: 450nm)
- 37°C incubator
- Automated plate washer
- Precision single and multi-channel pipette and disposable tips
- Clean tubes and Eppendorf tubes
- Deionized or distilled water

7. Washing

Manual:

- Discard the solution in the plate without touching the side walls. Clap the plate on absorbent filter papers or other absorbent material. Fill each well completely with 350ul wash buffer and soak for 1 to 2 minutes, then aspirate contents from the plate, and clap the plate on absorbent filter papers or other absorbent material.

Automatic:

- Aspirate all wells, and then wash plate with 350ul wash buffer. After the final wash, invert plate, and clap the plate on absorbent filter papers or other absorbent material. It is recommended that the washer shall be set for soaking 1 minute. (**Note:** set the height of the needles; be sure the fluid can be sipped up completely)

8. Sample Collection and Storage (universal)

- **Serum:** Place whole blood sample at room temperature for 2 hours or put it at 2-8°C overnight and centrifugation for 20 minutes at approximately 1000×g. Collect the supernatant and carry out the assay immediately. Blood collection tubes should be disposable, non-pyrogenic, and non-endotoxin.
- **Plasma:** Collect plasma using (EDTA-Na₂ or heparin as an anticoagulant. Centrifuge samples for 15 minutes at 1000×g at 2-8°C within 30 minutes of collection. Collect the supernatant and carry out the assay immediately. Avoid hemolysis, high cholesterol samples.
- **Saliva& Nasal fluid:** Centrifuge samples for 20 minutes at 10000×g at 2-8°C. Collect supernatant and carry out the assay immediately.

Note: Samples to be used within 5 days can be stored at 2-8°C, besides that, samples must be stored at -20°C (assay \leq 1 month) or -80°C (assay \leq 2 months) to avoid loss of bioactivity and contamination. Avoid multiple freeze-thaw cycles. The hemolytic samples are not suitable for this assay.

9. Reagent Preparation and Storage

Bring all reagents and samples to room temperature for 20 minutes before use.

- **Wash Buffer:**

If crystals have formed in the concentrate, you can warm it with 40°C water bath (Heating temperature should not exceed 50°C) and mix it gently until the crystals have completely been dissolved. The solution should be cooled to room temperature before use.

Dilute 30ml Concentrated Wash Buffer into 750ml Wash Buffer with deionized or distilled water. Put unused solution back at 2-8°C.

- **Preparation of Biotin-labeled Antigen Working Solution:**

Prepare it within 30min before experiment.

1. Calculate required total volume of the working solution: $50\mu\text{l} / \text{well} \times \text{quantity of wells}$. (Allow 55-60 μl more than the total volume.)
2. Dilute the SABC with SABC Dilution Buffer at 1:100 and mix them thoroughly. (i.e. Add 1 μl of SABC into 99 μl of SABC Dilution Buffer.)

- **Preparation of Biotin-labeled Antigen Working Solution:**

Prepare it within 30min before experiment.

1. Calculate required total volume of the working solution: $50\mu\text{l}/\text{well} \times \text{quantity of wells}$. (Allow 55-60 μl more than the total volume.)
2. Dilute the SABC with SABC Dilution Buffer at 1:100 and mix them thoroughly. (i.e. Add 1 μl of SABC into 99 μl of SABC Dilution Buffer.)

10. Assay Procedure

When diluting samples and reagents, they must be mixed completely and evenly. Before adding TMB into wells, equilibrate TMB Substrate for 30 minutes at 37°C. It is recommended to plot a standard curve for each test.

1. Bring all reagents to room temperature before use.
2. Label the sample wells, 2 Negative Controls, 2 Positive Controls and 1 blank well. Wash plate 2 times before adding sample and control (blank) wells!
3. Add 40 μ L sample dilution buffer to each sample well.
Add 50 μ L sample dilution buffer to blank well.
4. Add 10 μ L sample to each sample well.
Add 50 μ L Negative Controls and Positive Controls to set Controls well and gently tap the plate to ensure thorough mixing. Seal the plate with a cover and incubate at 37°C for 30 min.
5. Remove the cover, and wash plate 2 times with Wash buffer and each time let the wash buffer stay in the wells for 1-2 min.
6. Add 50 μ L Biotin-labeled Antigen to each well. Seal the plate with a cover and incubate at 37°C for 30 min.
7. Remove the cover, and wash plate 3 times with Wash buffer and each time let the wash buffer stay in the wells for 1-2 min.
8. Add 50 μ L of SABC Working Solution into each well, cover the plate and incubate at 37°C for 30 minutes.
9. Remove the cover, and wash plate 5 times with Wash buffer and each time let the wash buffer stay in the wells for 1-2 min.
10. Add 50 μ L of TMB substrate into each well. Gently tap the plate to ensure thorough mixing. Cover the plate and incubate at 37°C in dark within 10-15 min. And the shades of obvious blue can be seen in the Positive Controls. Blank well wells show no obvious color.
11. Add 50 μ L of Stop solution into each well and mix thoroughly. The color changes into yellow immediately.

12. Read the O.D. absorbance at 450 nm in a microplate reader immediately after adding the stop solution (Use the blank well to set zero).

11. Data Analysis

Calculation of Results

$$\text{Cutoff Value} = \text{NCx} \times 2.1$$

NCx: Mean Absorbance of Negative Control (when $\text{NCx} < 0.05$, Calculate as 0.05).

PCx: Mean Absorbance of Positive Control

1. Sample with absorbance values $<$ Cutoff Value are considered negative.
Sample with absorbance value \geq Cutoff Value are considered positive.
2. $\text{PCx} \leq 0.5$, the test is regarded as invalid, should be tested again.

Sample test data (for reference only)

Samples came from rehabilitation clients (1-2 months after recovery) of mobile cabin hospital. The plasma samples were diluted 1:5. TMB Color development time was 15 minutes at 37°C. $\text{NCx} = 0.109$

Rehabilitation clients (OD450)				Healthy volunteers (OD450)			
1#	1.850	9#	1.876	1#	0.101	9#	0.112
3#	1.371	10#	1.556	2#	0.094	10#	0.095
3#	1.962	11#	2.010	3#	0.077	11#	0.088
4#	1.984	12#	1.458	4#	0.09	12#	0.099
5#	1.985	13#	1.957	5#	0.143	13#	0.134
6#	1.885	14#	2.144	6#	0.128	14#	0.120
7#	2.046	15#	2.202	7#	0.112	15#	0.122
8#	1.295	16#	1.891	8#	0.096	blank	0.139

Who we are

Attogene is a biotechnology company located in Austin, Texas. Our focus is to enhance health and wellness by offering and developing customer focused Life Science Products domestically and internationally.

Our mission is to:

- Enhance detection technologies
- Enable rapid responses
- Enable impactful research discoveries

Contact Us

3913 Todd Lane, Suite 310
Austin, TX 78744

Phone: 512- 333-1330

Email: sales@attogene.com

Web: www.attogene.com

EL2045-01.V2_20200714