



Certificate of Conformity

This certifies that the following product,


Microcystin ELISA Kit

Catalog No. EL2024-01 Lot No. ATT4236L
Release Date: 8/24/2021 Expiration Date: Aug/2022

Complies with Attogene's quality control standards and conforms to product release specifications.

Quality Control Assay Data

Calibrator Concentration (ppb)	Average Optical Density (nm)	Coefficient of Variation	B/B ₀	Microcystin Concentration (ppb)
0.00	2.13	1.74	N/A	N/A
0.10	1.70	3.34	80.1%	0.10
0.30	1.17	1.60	55.1%	0.30
0.80	0.76	0.32	35.7%	0.82
1.00	0.68	1.69	32.0%	1.07
2.00	0.54	0.98	25.3%	1.97



R ² Value
1.0