



Certificate of Conformity

This certifies that the following product,

Cylindrospermopsin ELISA Kit


Catalog No. EL2047-01 Lot No. ATT1174L
 Release Date: 07/01/2021 Expiration Date: Jun/2022

Complies with Attogene's quality control standards and conforms to product release specifications.

Quality Control Assay Data

	Rep. #1	Rep. #2	Rep. #3	Avg.	Std. Dev.	%RSD	%Bo	'ppb cylindro
0 ppb Calibrator	1.93	1.95	1.92	1.90	0.04	2.35		
	1.88	1.82	1.88					
0.05 ppb Calibrator	1.46	1.40	1.46	1.45	0.04	3.09	76%	0.05
	1.52	1.40	1.46					
0.3 ppb Calibrator	0.64	0.63	0.65	0.64	0.01	1.27	34%	0.30
0.75 ppb Calibrator	0.32	0.32	0.32	0.32	0.00	0.80	17%	0.77
2.0 ppb Calibrator	0.16	0.15	0.16	0.16	0.00	3.12	8%	1.95
0.25 ppb Control (±0.05ppb)	0.75	0.70	0.71	0.72	0.02	3.32	38%	0.25

Negative control absorbance must be greater than 1.2 OD and less than 3.0 OD at release date.



Attogene certifies that this lot has passed Quality Control assays and conforms to product release specifications.