



## Certificate of Conformity

This certifies that the following product,


### Brevetoxin Plate Kit

Catalog No. EL2023-01    Lot No. ATT1310L  
Release Date: 11/9/2021    Expiration Date: JUN/2022

Complies with Attogene's quality control standards and conforms to product release specifications.

### Quality Control Assay Data

Calibrator Concentration (ppb)	Average Optical Density (nm)	Coefficient of Variation	B/B <sub>0</sub>	Brevetoxin Concentration (ppb)
0.0	1.64	0.74	N/A	N/A
0.1	1.28	2.25	78%	0.10
0.25	0.97	1.20	59%	0.25
0.5	0.69	1.50	42%	0.52
1.0	0.50	2.88	30%	0.96
2.5	0.31	3.16	19%	2.60



R <sup>2</sup> Value
1.0