



Rapid Test Components Plastic Cassettes


| | |
|------------------------|--|
| Product | Plastic Cassettes for Rapid Tests |
| Intended Use | Development/Manufacture of Lateral Flow Immunoassay Test Kit |
| Product Variant | CO2022 |

| Property | All Type |
|------------|----------|
| Appearance | White |
| Fitment | OK |

| Type | Description |
|--------|-----------------------|
| CO2022 | 5 mm plastic cassette |

Specifications

REPRESENTATIVE DIAGRAM OF PLASTIC CASSETTE - NOT TO SCALE



| Parameters of Plastic Cassettes | CO2022 Dimensions |
|--|-------------------|
| A. Sample Hole Diameter | 4.1 mm |
| B. Width of Strip | 5.0 mm |
| C. Length of Strip | 60 mm |
| D1. Distance from sample hole to Pressure point (sample pad side) | 10 mm |
| D2. Distance from sample hole to Pressure point (Absorbent pad side) | 35.5 mm |
| D3. Distance between sample hole and lower side of square hole | 10.7 mm |
| L1. Width of the device | 20.5 |
| L2. Length of device | 70.6 |
| Length of window | 16.8 |
| Width of window | 4.2 |
| NC Membrane length on the strip | 25 |

Storage Condition Store and handle in accordance with all current regulations and standards. Store in original packing at room temperature. Avoid direct sunlight. **Caution:** Material should be kept away from heat radiators.