



# Lateral Flow Cassette

Catalog No.	Description
CO0013	71×20mm/63×3mm Cassette

## Physical Properties

Property	All Type
Appearance	White
Fitment	OK

## Storage Condition

Store and handle in accordance with all current regulations and standards. Store in original packing at room temperature. Avoid direct sunlight.

**Caution:** Material should be kept away from heat radiators.

## Specifications

Parameters of Plastic Cassettes	Dimensions
Width of Strip	3 mm
Length of Strip	63 mm
Width of the device	20 mm
Length of device	71 mm
Length of window	13.8 mm
Width of window	2.7 mm
Sample Hole Diameter	2.8 mm



# Lateral Flow Cassette

